

# Methods

## Sequencing

### 1. Extraction of genome DNA

Total genome DNA from samples was extracted using CTAB/SDS method. DNA concentration and purity was monitored on 1% agarose gels. According to the concentration, DNA was diluted to 1 ug/ $\mu$ L using sterile water.

### 2. Amplicon Generation

16S rRNA/18S rRNA/ITS genes of distinct regions(16S V4/16S V3/16S V3-V4/16S V4-V5,18S V4/18S V9,ITS1/ITS2, Arc V4) were amplified used specific primer(e.g.16S V4:515F-806R,18S V4: 528F-706R, 18S V9:1380F-1510R,et.al) with the barcode. All PCR reactions were carried out with 15  $\mu$ L of Phusion® High-Fidelity PCR Master Mix (New England Biolabs); 0.2  $\mu$ M of forward and reverse primers, and about 10 ng template DNA. Thermal cycling consisted of initial denaturation at 98°C for 1 min, followed by 30 cycles of denaturation at 98°C for 10 s, annealing at 50°C for 30 s, and elongation at 72°C for 30 s. Finally 72°C for 5 min.

### 3. PCR Products quantification and qualification

Mix same volume of IX loading buffer (contained SYB green) with PCR products and operate electrophoresis on 2% agarose gel for detection. PCR products was mixed in equidensity ratios. Then, mixture PCR products was purified with Qiagen Gel Extraction Kit(Qiagen, Germany).

### 4. Library preparation and sequencing

Sequencing libraries were generated using TruSeq® DNA PCR-Free Sample Preparation Kit (Illumina, USA) following manufacturer's recommendations and index codes were added. The library quality was assessed on the Qubit®2.0 Fluorometer (Thermo Scientific) and Agilent Bioanalyzer 2100 system. At last, the library was sequenced on an Illumina NovaSeq platform and 250 bp paired-end reads were generated.

## Data analysis

### 1.Paired-end reads assembly and quality control

#### 1.1 Data split

Paired-end reads were assigned to samples based on their unique barcode and truncated by cutting off the barcode and primer sequence.

## 1.2 Sequence assembly

Paired-end reads were merged using FLASH (V1.2.7, <http://ccb.jhu.edu/software/FLASH/>)<sup>[1]</sup>, a very fast and accurate analysis tool, which was designed to merge paired-end reads when at least some of the reads overlap the read generated from the opposite end of the same DNA fragment, and the splicing sequences were called raw tags.

## 1.3 Data Filtration

Quality filtering on the raw tags were performed under specific filtering conditions to obtain the high-quality clean tag<sup>[2]</sup> according to the QIIME (V1.9.1, [http://qiime.org/scripts/split\\_libraries\\_fastq.html](http://qiime.org/scripts/split_libraries_fastq.html))<sup>[3]</sup> quality controlled process.

## 1.4 Chimera removal

The tags were compared with the reference database (Silva database, <https://www.arb-silva.de/>) using UCHIME algorithm (UCHIME Algorithm, [http://www.drive5.com/usearch/manual/uchime\\_algo.html](http://www.drive5.com/usearch/manual/uchime_algo.html))<sup>[4]</sup> to detect chimera sequences, and then the chimera sequences were removed<sup>[5]</sup>. Then the Effective Tags were finally obtained.

# Reference

- [1] Magoc T, Salzberg S L. FLASH: fast length adjustment of short reads to improve genome assemblies. *Bioinformatics* 27.21(2011):2957-2963.
- [2] Bokulich, Nicholas A., et al. Quality-filtering vastly improves diversity estimates from Illumina amplicon sequencing. *Nature methods* 10.1 (2013):57-59.
- [3] Caporaso, J. Gregory, et al. QIIME allows analysis of high-throughput community sequencing data. *Nature methods* 7.5 (2010):335-336.
- [4] Edgar, Robert C., et al. UCHIME improves sensitivity and speed of chimera detection. *Bioinformatics* 27.16 (2011): 2194-2200.
- [5] Haas, Brian J., et al. Chimeric 16S rRNA sequence formation and detection in Sanger and 454-pyrosequenced PCR amplicons. *Genome research* 21.3 (2011): 494-504.